

The PROGRESSOR demo and slide show presentation download contains a PowerPoint presentation which explains the theory and function of the PROGRESSOR software, and a demo version of the PROGRESSOR software itself.

Download the file progressor\_demo.zip from

[www.medisoft.co.uk/progressor/progressor\\_demo.zip](http://www.medisoft.co.uk/progressor/progressor_demo.zip)


and extract this file to a destination on your hard drive, or other device such as a USB memory source.

To extract the file using Windows XP, right-click on the file and select 'Extract All...'. You will be prompted to choose a destination location to extract the files to. If you do not have this option (e.g. Windows 2000) you will need to install WinZip or an equivalent unzipping program, in order to extract the files ([www.winzip.com](http://www.winzip.com))


Once you have extracted the zip file to your destination location, browse to that location, for



example, by double-clicking on your 'My Computer' icon and then navigating to your destination location. In the example below this is the D:\ drive. Double click on the file called 'Play' which you will subsequently see as highlighted in the following example:



'PowerPoint Viewer' will start. You may be asked to accept a licence agreement:




If so, simply click 'Accept' and the presentation will start.

You do not need to have PowerPoint installed on your computer. The viewer, presentation and PROGRESSOR Demo all run directly from disk, and no files are installed on your hard drive.

Once the presentation starts, you will have 4 options:

- Demo** Click to run PROGRESSOR as a demo version.
- Show** Click to view a presentation explaining PROGRESSOR (approx 10 mins).  
NB: the presentation has a commentary throughout, so make sure your volume / speakers are switched on if you computer has speakers.
- Buy** Click for an e-mail link to contact Medisoft about purchasing PROGRESSOR.
- Quit** Exit the demo and presentation.

When you click the 'Demo' button you may be asked the following question:



If so, click 'Yes' to run the demo. This may then be followed by another question:



If you see this message, click 'Run'.

When PROGRESSOR starts you will see a screen with 2 example patients

Patient database

Name	ID
PATIENT 1	1111
PATIENT 2	2222

Analyze

Combine

Data import

Update


Delete

Rename

Help

Auto update

Click on one of the patients and then click the 'Analyze' button. This brings up the progression analysis for the patient, e.g.



Consult the contents of on line help file, and watch the slide show presentation to learn about the full list of features for analysing visual field data in PROGRESSOR.