



Creating Surgeon Defaults

Always review data loaded by default to ensure it is correct.


1 Log In
Open MediSIGHT. Enter your Username and Password. Click **Login**.




2 Open System Settings
Click on **System Settings** at the top of the Home page.




3 Open Defaults
Click on **Defaults** at the top of the Settings page.




4 Open Drug Defaults
Click on **Drug Defaults - User**
Select the user.
Select the description.




5 Add Drug Details
Click on the **Add** button.
Complete the Add Drug page.
Click on **Save**.



6 Edit Drug Details
Use the Edit icon to make changes to a drug and **X** to delete it.
Click **Go Back** to return to the Defaults page.



7 Open Surgeon Defaults
Go to **Surgeon Defaults** and click **Create surgeon default**.
Select the procedure and click **OK**.



8 Enter and Save Surgeon Defaults
Complete each tab.
Click **Save & Close**.

