


Creating a Follow-up Visit


1 Log In
Open mediSIGHT. Enter Username, Password then click **Login**.




2 Find Patient
Click on the **Clinic List** step or type the Hospital/NHS Number in **Quick** search. Select the search result to open the record.



3 Create Follow-up
Check the patient details at the top before you start.



Click on **Create new +**
Select **Cataract follow-up**



4 Enter Data
Click on a section to expand it. Click **Done** as you enter data to collapse that section. You must complete all *fields.

Pre-operative assessment

Pathway * Macular Treatment Location * Eye Clinic (WH)

Date * 18-May-2018 Responsible


Seen by * Ms Sarah LEWIS

In clinic By telephone


Done

Is data entered for the correct eye?


5 Save and Print
Use the **Finalise** button to save and confirm the data. Print-outs will be sent to your default printer.



6 Check Finalised Data



7 Edit Data
Find the event and use the crayon if you need to edit the data.



8 Log Out
Use your profile drop-down menu to **Log out** when finished.

